



Cochlear™ Wireless Guide

Use this guide to help select Cochlear™ Wireless Accessories.



Cochlear Wireless Phone Clip

Select this accessory to:

- Use hands-free calling directly to the sound processor
- Use a wireless headset to make phone calls on Bluetooth® enabled phones
- Hear Bluetooth enabled navigation system directly in the sound processor
- Stream music from a Bluetooth enabled device

The Phone Clip clips onto clothing allowing the built-in microphone to pick up the user's voice, while the phone can be several meters away.



Cochlear Wireless Mini Microphone

Select this accessory if you want:

- To hear better in noisy listening situations:
 - In a meeting, riding in a car, dining at a restaurant, or participating in a class at the gym
- A portable streamer for watching TV at a friend's house
- To stream music from an iPod® or music streamer

The Mini Microphone can be clipped onto clothing and will transmit speech wirelessly over a distance up to 6 metres. Multiple users can be connected to the same Mini Microphone.



Cochlear Wireless TV Streamer

Select this accessory to:

- Enjoy TV viewing without turning up the volume and disturbing others
- Listen to TV, stereo, or other audio source
- Move around freely while listening to the TV

The TV Streamer allows multiple users to connect to the same device, bringing together entire families of users. 2.4 GHz technology means there is little chance of delay between what you see and what you hear.

NOTE: Cochlear has looked at how we use our processors and has developed three wireless accessories to help us. They are only compatible with the N6 processors and Baha. For a Cochlear Implant they do need to be paired to your processor remote and must be purchased through your Cochlear Implant Clinic.



Hear now. And always.

